

Testing Solutions for pMDI Drug Development



Pressurised metered dose inhalers (pMDIs) are a class of combination drug products dependent on the optimisation of a formulation, device design and human usage to deliver an accurate, reproducible dose. In the case of pMDI suspension products, the shaking profile is crucial for accurate dose delivery. Lack of appropriate shaking can deliver a high amount of the drug in the early doses followed by very little or no drug towards the end of the product's life. Other important influencing factors on pMDI performance include force-to-actuate and hold time to ensure that the metering valve is open long enough to deliver a complete dose. A minimal shake-to-fire delay is another crucial factor for maintaining the uniformity of the delivered dose (Figure 1).

COMPLETE SUITE OF pMDI TESTING SOLUTIONS

Proveris by Design

Based on the quality by design (QbD) approach endorsed by the US FDA, Proveris by Design first provides the basis of experiments to identify the range of human usage parameters of a pMDI according to current regulatory guidelines for *in vitro* testing. Next, it offers strategies to test the range of design spaces and

“The insight gained into product performance can be used long after product approval to investigate and resolve any out-of-specification or out-of-trend disruptions to manufacturing.”




Figure 1: Example of dose content and shot weight for a suspension pMDI product with improper and proper shaking.

gain an understanding of the control space: which of these factors have the most impact on product performance? Lastly, it identifies the target operating space and control for key influencing factors in a

tight range during the *in vitro* testing. One of the deliverables of the approach is a sensitivity map based on the design of experiments which shows how much influence certain factors identified

Typical Influencing Factors for pMDI Products		Screening Experiments Category
Formulation	Morphology of API (size and shape) Formulation form (suspension/solution) Propellant Excipients	Performance testing
Device	Metering valve Actuator (sump geometry) Assembly process	Component characterisation
Human Usage	Shaking and shake-to-fire delay Hold time Product actuation force Actuation velocity	Functional testing

Table 1: Typical influencing factors for pMDI products and Proveris by Design screening experiments.

in the design space have on the overall performance of the product. In the R&D phase, clients engage Proveris Laboratories to perform tests using this method to produce consistent, high-quality data efficiently during product development and for preparation of regulatory submissions. The results show an overall reduction in approval timeline and minimised queries from regulatory agencies. Moreover, the insight gained into product performance can be used long after product approval to investigate and resolve any out-of-specification or out-of-trend disruptions to manufacturing (Table 1).

Ergo Studies

As a contract test service, Proveris Laboratories offers human-realistic actuation studies performed using proprietary Ergo technology to quantify accurately how trained testers in the targeted population actuate the product. Client device candidates, and/or reference products in the case of generic drug development projects, are evaluated to obtain key actuation parameters, such as stroke length, velocity, acceleration and hold time. This data is fully transferrable to Proveris Vereo automated actuators to enable reproducible actuations in a human-realistic way.

Vereo SFMDx Actuators

As with any pharmaceutical product, it is important to perform *in vitro* tests with repeatable, reproducible and robust methods for a smoother product development process. Vereo actuators provide fully controllable and repeatable shaking and actuation for commonly required regulatory tests, controlling up to six critical actuation parameters and four additional parameters related to shaking for pMDI devices. The flexible SFMDx actuator fits seamlessly into multiple testing workflows with identical actuation parameters implemented across all tests, ensuring data integrity and accuracy (Figure 2).

In vitro tests performed using Vereo SFMDx actuators (Figure 3) include:

- Delivered dose uniformity (DDU)
- Single actuation content uniformity (SAC)
- Valve/pump delivery (shot weight)
- Aerodynamic particle size distribution (APSD)
- Droplet size distribution (DSD)
- Spray pattern and plume geometry

Force - Time - Position Graph




Figure 2: Proveris instrument users can create and evaluate products using force/position versus time plots generated by the system software.



Figure 3: Proveris precision instruments and test workflows incorporate Vereo actuator technology for testing of pMDI products.

- Priming and repriming
- Product wasting for through-life testing

testing, using a powerful common software platform consistent with full regulatory compliance (21 CFR Part 11).

New Kinaero Cx pMDI Collection System

The newly introduced Kinaero Cx system easily integrates into DDU and APSD test workflows to automate parameters including shaking, actuation and dose/sample collection in a reliable, repeatable manner using proven Vereo actuator technology. To automate the whole workflow, users can run the Kinaero Cx alongside the Kinaero High-Throughput pMDI Fire-Down System for through-life

Proveris Scientific Corporation

Two Cabot Road
Hudson
MA 01749
United States

T: +1 508 460 8822
E: contactus@proveris.com

www.proveris.com

NEW!

kinaero cx™

pMDI Collection System



ACHIEVE ACCURATE, REPRODUCIBLE RESULTS FOR DDU & APSD TESTING

MAXIMIZE PMDI TESTING EFFICIENCY

Introducing a new addition to the Proveris product family!

Building on more than 25 years of innovation, expertise, and specialization in inhalation product analysis and performance testing, the **Kinaero Cx pMDI Collection System** automates shaking, actuation, and dose collection for delivered dose uniformity (DDU) and aerodynamic particle size distribution (APSD) testing.

- Operation and data management software with database storage and retrieval ensures regulatory compliance (FDA 21 CFR Part 11)
- Precise control of shake-to-fire timing, actuation parameters, and flow rate
- Gold standard Vereo® Automated Actuator technology
- Flexibility to run DDU and APSD measurements on one platform
- System works with most pMDI devices

See the Kinaero Cx System in action
at RDD 2022 in Orlando, Florida

Email contactus@proveris.com for more
information or stop by our table

LEARN MORE:

www.proveris.com
contactus@proveris.com
+1 508 460-8822

